

CPL-2010 COLD PLATE



- The high performing unit is cooling down to -30°C in a very short time.
- The electronic temperature control with a range from -5°C to $+30^{\circ}\text{C}$ is continuously adjustable.
- The following accessories, available at moderate cost, add versatility to the application of this high performance quick cool unit.
- The temperature is adjustable and visible by digital control.
- The cold plate is made of polished stainless steel, with heated brims to prevent condensation
- There is also acrylic lid to prevent condensation.
- The system can be use with dispenser unit.

Technical Specifications

Dimension(WxDxH) : 652x380x425 mm

Weight : 20kg

Working condition : Room temperature, it can be work until %60 relative humidity.

Power supply : 20V \pm % 10 AC, 50 Hz



DES San. Sit. 1.Cd. D3 Bl. No:27 34776
Yukari Dudullu, Istanbul - TURKEY

Phone : +90 216 365 40 99

Fax : +90 216 365 40 89

info@baygenlab.com
www.baygenlab.com