

SF-44 GROSSING WORKSTATION



- All stainless steel construction (1.4301-ANSI 304), working area (ANSI 316)
- Working area can be perforated
- All functions working by photo sensor and manual
- Backdraft-Downdraft ventilation removes formaldehyde fumes safely
- Formalin waste tank
- Advanced electrical system



■ Standard Accessories

- Paper towel holder
- Removable measuring rule
- Magnetic instrument holder
- Polyethylene dissecting board
- Glove box holder
- Heavy-Duty disposal unit
- Fluorescent/Led light
- Formalin waste tank
- Magnifier
- Formalin faucet
- Hot/Cold water faucet




baygen

Dimensions W/D/H	:	150/75/200 cm
Exhaust Air capacity	:	1600 m3/h
Construction Material	:	1.4301-ANSI- 304-306
Working area	:	142*65 cm
Cold/Hot water connection	:	½ "
Waste water connection	:	2"
Model no	:	SF 44
Custom manufacture	:	Applicable

Sensor & Button Group

All functions working by photo sensor and also manual

- Cold water
- Emergency/stop
- Hot water
- Fan
- Surface washing
- Disposal unit
- Leakage relay
- Light



■ Optional Accessories

- Automated formalin dispensing system
- Cpu stand
- Cassette holder
- Computer monitor and keyboard stand
- Splash shield
- Camera mount facilities
- Video camera facilities
- Dictation equipment stand
- Glove box holder
- Dictaphone
- Document holder



Automated formalin dispensing system



Dictaphone



Computer monitor and keyboard stand

DES Sanayi Sitesi 1. Cadde
D-3 Blok No:27 34776
Dudullu Istanbul

Tel : +90 216 365 40 99

Fax : +90 216 365 40 89

info@baygenlab.com

www.baygenlab.com/en

baygen